

CENTRIFUGE INSTRUCTIONS

Everyone must be trained on the centrifuge before they can use it.

The centrifuge is NOT available after hours.

1. DO NOT use centrifuge if you are unfamiliar with it. Speak with the contact person first for a demonstration.
2. Sign out centrifuge every time you use it.
3. Return rotor to cupboard below centrifuge. DO NOT place them on lab benches.
4. Clean out any spills in rotors immediately. Remember to consider what you are working with when disposing of spilled materials.
5. Do NOT force screw into rotor and shaft. If you are having difficulty, **STOP** immediately, unscrew, lubricate and try again. Ask for help if you are having difficulties.
6. Lubricate screw periodically to avoid any mis-threading.

CONTACT: Diane Kinney (Rm. 202, Ext. 2859)