

## **IMAGING COMPUTER**

1. **DO NOT** wear gloves when using the iMac.
2. The priority of the imaging computer is for gel analysis and imaging. If you are using the computer for other purposes please be respectful to the fact that you may be interrupted.
3. The printer that should be used for gel images is the Brother HL-5070N (“PCR lab printer”) located in Rm 231. Please **DO NOT** change the printer default.
4. When printing gels:
  - a. Invert image
  - b. Decrease size of image
5. Save all files to the appropriate folder. Files remaining on the desk top will be deleted.
6. Do NOT install any software or download any programs without receiving permission.

**Contact: Melanie Lewis Ivey (Rm 231, ext 2704)**